



QSM Quantum Scanning Microscope

NV-MAGNETOMETER: FASTER. VERSATILE. EASY TO USE.

Our NV Magnetometer and its Applications

Meet Our Turn-key Scanning NV Microscope: The QSM

The QSM unites user-friendliness with flexibility and high performance. From large area megapixel images to high resolution, high sensitivity maps, the QSM is the right tool for the job. Delivering perfomance close to the shot noise limit, its flexible platform makes it a versatile tool for studying magnetic phenomena at the nanometer scale, from stray fields and vortices over domain boundaries to current densities.

Flexibility at the Core of Our Design



Measurement for Scanning Modes

NV-CONTINUOUS WAVE (CW)

Oualitative: PL. Iso-B Quantitative: ODMR



NV-PULSED

Pulsed ODMR, Ramsey, CPMG, arbitrary user-defined



MOKE

Magneto-optical Kerr Effect



ADDITIONAL

AFM topography with height control (two-pass mode and fixed plane)

Widefield & confocal

Multiple experiments per pixel possible (example: ODMR, Rabi, CPMG)

Working Principle

We use a single defect in diamond, a Nitrogen Vacancy (NV) center, as an atomic size quantum sensor. By scanning it in close proximity over a sample surface, we can quantitatively map magnetic fields and more with high



Applications at a Glance

NV SENSOR FOR

PURPOSE

Magnetic Fields

Currents

Electrical Fields



Failure Analysis Materials Characterization Nano NMR Radical Sensing **RF SOIC R&D** Surface Characterization Thermometry Transistor R&D Waveguide Characterization



sensitivity and better than 40 nm resolution. All of this works in ambient conditions, in the desktop sized QSM.

Removing two adjacent carbon atoms from the diamond crystal lattice and placing a nitrogen atom at one of the sites creates an NV center. This "artificial atom" has several unique properties:

- It is a room temperature qubit •
- The levels can be readout via an optical signal
- The gyromagnetic ratio of 28 GHz/T makes it very sensitive to magnetic fields
- The measured magnetic fields are absolute

MATERIALS / DEVICES

2D Materials Antiferromagnets Hydrodynamic Tr. lon Traps Multiferroics **Optoelectronics Phase Transitions** Skyrmions **Spintronics** Surface Potentials **Topological Materials**

- Use Cases ·

ANTIFERROMAGNETS Domains in CuMnAs



Resolve the challenging, weakly varying magnetic signals.

Correlate magnetic field and magnetization

Quickly scan large sample areas to identify

for a better understanding of the material.

features of interest and use the additional data

MOKE

SKYRMIONS

Vortices with applied magnetic field



Tune the magnetic field to induce vortex nucleation.

TRACING CURRENTS



MULTIFERROICS

Cycloids in strained BiFeO

Study spin cycloid domains with

high sensitivity and resolution

Image currents in a running device via the magnetic fields they generate to reveal defects and hot spots or verify the design.

QUICKSCAN MODE Stripe features in Hematite (Fe₂O₃)



Image large areas while retaining high pixel density in a short amount of time with Quickscan. Resolve both short- and long range features.

RELAXOMETRY

Different spin decay rates due to surface-magnons coupling in YIG



The QSM can achieve Rabi rates exceeding 40 MHz.



Extend the QSM with multiple microwave sources, for example to perform spin selective readout.

SPECTROSCOPY

Ferromagnetic Resonance in YIG



Measure at high magnetic fields with our magnet options and broad microwave band.

Correlate topography with currents and more





DOMAIN WALLS Extract domain wall chirality and width



Characterize domain walls directly from the data, which is the projection of the sample magnetic field onto the defined NV axis.

NV ORIENTATIONS

In-plane, out-of-plane and standard orientation



We offer different NV orientations for optimum projection of the stray field or measuring at high bias.



System Overview -



OPTICS HEAD

Scanning confocal optics for NV excitation and readout, scanning MOKE and widefield imaging

AFM STAGE

AFM head with manual 3D stage for microwave excitation for NV scanning probes

MAGNET STAGE

Mechanical magnet or electro vectormagnet for biasing the NV and/or sample

SUPPORT RACK Quantum Control Console,

AFM controller, PC, UPS

MONITOR

Display real-time temperature and humidity data over up to 24 hours

ISOLATING ENCLOSURE

Temperature stabilized and acoustically isolated enclosure, active vibration isolation



QSM Highlights -

MULTI-COLOUR EXCITATION

The QSM optical setup and its control electronics can handle multiple wavelengths, as well as multi-colour excitation pulses and pulse sequences.



MAGNET OPTIONS

Mechanical magnet for bias fields exceeding 200 mT or vector electromagnet for fast sweeps.



Flexibility

MOKE MEASUREMENT SETUP

MOKE provides additional data in the form of sample surface magnetization and allows quicker identification on a larger size scale.

User-Friendly

UNIFIED USER INTERFACE

The QS3 software interface gives you full control of the QSM microscope, from the AFM to the NV part. It features unique measurement protocols such as Quickscan ODMR, is open source, and allows scripting within the same GUI.

BRIGHT FIELD IMAGING WITH HIGH RESOLUTION AND LARGE FIELD OF VIEW (FOV)

By switching between these two independent systems, you can easily identify larger and smaller landmarks on your chip. Approaching the surface is safe and easy.

> Green arrow: A 50x high resolution system with a FOV of 150 µm Blue arrow: A low resolution system with a FOV of ~2 mm

OPTIMIZED OPTICS

Our custom optical design significantly boosts the optical performance. The NV emission overlaps the visible and near IR region, which means standard components are not optimal. Each element in the QSM is optimized for

the NV band, even the objective.



BROADBAND MW EXCITATION FOR A LARGE RANGE OF MAGNETIC FIELD BIASES

5 GHz of microwave bandwidth means that the QSM can measure at any magnetic bias field up to 110 mT. Both magnet options allow full vector control of the bias field, giving you maximum Performance measurement flexibility.

QUICKSCAN

ODMR at pixel rates up to 200 Hz means 100 x 100 pixels scans in less than one minute and Megapixel scans in less than 2 hours. Optimized software algorithms decrease software overhead close to zero.



Powerful Software Suite

Background Information -

QS3 Software - Scanning NV as simple as AFM

What you see with your eyes is just the tip of the iceberg. The QS3 software is built around the experience and the needs of our customers combining functionality, speed and usability. The intuitive software interface jump starts your entry into the world of magnetic and current imaging. A unified, turnkey solution, the software seam-lessly integrates all the QSM components from the AFM to the magnetometer and the MOKE.



Quickscan - Megapixel Images in Less Than Two Hours

Measure faster with Quickscan! A unique algorithm practially eliminates the software overhead, so that the unmatched performance of our probes and optimized optical design of the QSM directly translate to faster measurements. At pixel rates of up to 200 Hz, 100 x 100 pixel scans don't take hours anymore but less than a minute. You can correlate short- and long-range features by taking megapixel scans in a few hours instead of weeks. Quickscan measures full frequency spectra at each pixel, providing the same wealth of information as conventional cw-ODMR.

Do you have a sample with highly varying field? Quickscan can also frequency track, keeping a small frequency window dynamically centered on the ODMR resonance. The small window preserves the sensitivity, which would be lost with the large frequency window necessary without tracking, keeping your scan quick!

Scripting and Pulsed Measurements

Being at the forefront of research means new and often unconventional experiments. A dedicated pulsed program interface in QS3 allows you to easily modify existing pulse sequences or even add new ones. Script custom pulse experiments controlling both the laser and the microwave source and directly use them to perform scanning measurements without additional programming.

A dedicated widget gives you easy access to software resources minimizing programming overhead. QS3 automatically plots and saves the data, so that you can focus on designing your experiments. Have a look at the spectroscopy measurement of page 5 to see what the script widget can do!

Sensitivity

Magnetic sensitivity strongly depends on the protocoll used. To start, let's look at cw-ODMR, measuring at a single point frequency point. The shot noise limited sensitivity is then given by: η_{cw}

$$\eta_{cw} = \frac{4}{3\sqrt{3}} \frac{h}{g_e \mu_B} \frac{\Delta v}{C_{cw} \sqrt{R}}$$

with linewidth Δv , contrast C_{cw} , count rate R.

In cw-ODMR, both contrast and count rate can be increased at the expense of linewidth. Optimal sensitivity is typically achieved at low linewidths rather than high contrast or count rates. **What values can you actually expect?** Our top of the line sensors have a contrast of better than 25% at a linewidth of 10 MHz and a countrate of up to 1 Mcts/s at maximum contrast, leading to a theoretical sensitivity of 1 μ T/ (Hz^{1/2}) in cw ODMR.

850

750

R [kcts/s]

nt rate

For real world instruments, soft- and hardware overhead as well as setting can significantly decrease the actual sensitivity. With Quickscan and tracking, the QSM can operate close to the shotnoise limit.

Measurement Modes Explained

PHOTOLUMINESCENCE

Magnetic fields modify the fluorescence rate of the NV, indicating where the field is stronger and where it is weaker. Since it only measures the fluorescence, it is a very fast, but only qualititave measurement.

ISO-B

Measures the ratio of two microwave frequency points centered around a specific frequency and about the line width apart to create contours of equal field. A fast measurement that is quantitative for small field changes.

CW-ODMR

Measures the full resonance curve by continually applying the laser and microwave field to map the fluorescence response. Fully quantitative with additional information in the contrast, countrate and linewidth.

TRACKING

Strong variations in the magnetic field shift the resonance significantly. This either requires a large frequency window for ODMR, which decreases signal-to-noise, or, to avoid this, a dynamically moving, small window.

PULSED ODMR

At each frequency point the NV is initalized with a laser pulse, followed by a microwave pi pulse and laser readout pulse. This avoids power broadening, the more narrow linewidth and higher contrast increase the sensitivity to better than 0.5 μ T/ (Hz^{1/2}).

GRADIOMETRY

The tuning fork motion means that the NV oscillates relative to the sample. This turns gradients into AC magnetic fields, which are detected with a low noise spin echo meausurement, increasing the sensitivity down to the 100 nT/ $(Hz^{1/2})$ range.



System Extensions -

Analog and Digital Extension

FASTER DATA. DIFFERENT DATA.

Adding the hard- and software for pulsed measurements significantly extends the capabilities of the QSM:

- Quickscan for up to 200 pixel/s scanning rates
- Pulsed ODMR, Rabi, Ramsey, T1, relaxometry, CPMG and more standard protocols integrated into the UI
- Pulse sequence generator for arbitrary, user-defined pulse sequences
- All pulsed protocols can be run stand-alone or as maps
- Protocols can be "stacked" for mapping. Example: ODMR/ Rabi/Ramsey to stay on resonace.

Pulsed measurements increase sensitivity and extend the measurements from DC magnetic fields to AC magnetic fields and currents, magnetic noise, relaxometry, gradiometry and much more. The additional up to 4 analog and 17 digital channels can be further used to apply stimuli or trigger additional instrumentation or hardware.

Magnetic Field Control

at bias fields beyond 100 mT.

OPTION 1: MECHANICAL MAGNET



OPTION 2: VECTOR ELECTROMAGNET

The low noise, temperature stabilized design offers >80 mT in arbitrary directions with peak fields for specific directions >200 mT.

The vector magnet enables easy field sweeps and ramps for magnetizing samples, as well as optimization algorithms for maximum NV performance.



Custom PCBs

OPTIMIZED TO YOUR NEEDS

Whether you need current pulses to move a domain wall or microwave excitation to excite a ferromagnetic resonance, applying signals to a sample is often required. We offer both standardized as well as custom designed ceramic PCBs with DC and/or RF lines, including connectors.

Contact us for features such as multiple microwave or laser sources, vacuum, temperature range, additional optical beam paths, rotating sample holder or inert atmosphere housings.

Specifications -

SPM Scanner and Coarse Stage

SAMPLE SCANNER	xyz scanner range: 85 μ System noise z: <50 pm
COARSE STAGE	xyz range 6 mm x 6 mm Stitching with scanner n
SPM MODES	Tuning fork based AFM Quickscan, pulsed proto
SAMPLE SIZE	25 mm diameter, (up to
PROBE-SAMPLE DRIFT RATE	<2 nm/h over 50 h in ter

Imaging Modes

NV SCANNING	Magnetic spatial resolut Magnetic sensitivity*: 0.
AFM TOPOGRAPHY	Tip dependent: NV scan Akiyama probes (non N
CONFOCAL FLUORESCENCE	Circular scanning range
ADDITIONAL	MOKE, bright field imagi

* Probe and protocol dependent: shot noise limit for DC sensitivity in cw-ODMR : better than 1.5 µT/(Hz^{1/2}), for pulsed ODMR better than 500 nT/(Hz^{1/2})

Confocal Microscope Performance

MICROSCOPE OBJECTIVE	50x magnification, NA=0 NV-band optimized and
OPTICAL SCANNING RANGE	Circular aperture with 15

Ouantum Control Console

MICROWAVE SOURCE	0.6-6 GHz, up to 36 dBm (5 W), cu Analog and digital extension: AWC Custom MW circuitry available on
LIGHT SOURCE & DETECTION	515 nm laser source (cw & pulsed SPAD module with high quantum

Vibration and Temperature Isolation

VIBRATION ISOLATION	Active isolation 0.7-300 Transmissibility <0.01(~
TEMPERATURE STABILIZATON	<0.1°C temperature drift



m x 85 µm x 15 µm (closed-loop, 0.15 nm resolution) RMS (dependent on environment)

x 15 mm (100 nm resolution) novement possible

mode, NV modes: Photoluminescence, Iso-B, cw-ODMR, cols

50 mm x 50 mm with custom holder)

mperature stabilized housing

tion: theoretical <10 nm; typical 25-70 nm 5-3 μ T/(Hz^{1/2})

ning probes \sim 1 nm (z), \sim 150 nm (xy); V) <1 nm (z), <30 nm (xy)

: 150 μm diameter; <1 μm imaging resolution @515 nm

ing (150 μm and 2 mm FOV)

).75, WD=5 mm, FOV=150 μm diameter, compatible with 515 nm, 575 nm and 637 nm

50 µm diameter FOV

stom range and output power available on request G with 1.25 GS/s, 400 MHz bandwidth request

l control), >10 mW, <1.5 ns rise/fall time efficiency

Hz, passive isolation above 300 Hz 40 dB) above 10 Hz



SHAPING A NEW FUTURE FOR NANOSCALE IMAGING

Established in 2018, QZabre AG is a pioneer in NV technology. As an ETH Spinoff from the Spin Physics Group, our vision is to make scanning NV imaging as easy to use as AFM today.



65

Regina-Kägi-Strasse 11 CH-8050 Zurich

www.QZabre.com info@QZabre.com

